

CPM Instructions

The CPM was prescribed by your physician to assist you in your rehabilitation. If you have any questions about your device please call your JACE Systems Sales Representative. If you have questions regarding your CPM therapy please call your therapist.

For your safety before starting your CPM therapy:

1. Check knobs to be sure they are tight.
2. Check your ranges and speed to be sure they are properly set. Remember, speed is set to your comfort.
3. Check to be sure the device is secure. Not “sliding or migrating”.
4. Always begin and end your CPM therapy at a low and comfortable range.
5. Begin by activating the patient switch.
6. Keep the patient switch in your hand so you can stop the device if you feel any discomfort.

JACE Hand and Wrist - TO CHANGE PROGRAM

1. Stop motion in comfortable position.
2. Press 'PR' arrow twice.
3. Press 'UP' or 'DN' arrow to desired program.
4. Press 'NX' arrow twice. Push GO.

TO CHANGE RANGE OF MOTION

1. Stop motion in comfortable position.
2. Press 'PR' arrow once.
3. Press 'UP' arrow once.
4. Press 'NX' arrow once.
5. Flex limit will display. Press 'UP' or 'DN' arrow to desired limit.
6. Press 'NX' arrow. Extn Limit will display. Press 'UP' or 'DN' arrow to desired limit.
7. Press 'PR' arrow (4 times) until screen displays 'Push GO To Start'.
8. Push 'GO'.

NOTE: If the display ever reads 'SET ORTHO UP', press and hold the up arrow until it reads 'Push GO to Start'.

JACE Ankle/Knee/Elbow

1. Turn the “Power” switch on
2. Begin and end your therapy at a comfortable range
3. Position your leg, foot or elbow as instructed. Positioning is very important and should be consistent for every session.
4. Attach any straps you may be using around your foot and ankle
5. Begin by activating your patient switch, keeping it nearby for easy access.

JACE Shoulder

1. Check black knobs on adjustment bars for position as well as tightness
2. Check to make sure device attachment to chair is secure and at proper height
3. Make sure correct motors are plugged in
4. Turn the “Power” switch on
5. Begin and end your CPM therapy with the arm in its lowest position
6. Position your shoulders so that they are level and relaxed
7. Sit with your back straight and flush up against the chair
8. Begin by activating your patient switch, keeping it nearby for easy access.

JACE Systems Corporate Customer Service Office hours: 9AM-6PM EST - 800-800-4276.

Thank you for using JACE Systems!