



JACE W550 Wrist CPM Instructions

The JACE W550 was prescribed by your physician to assist you in your rehabilitation. If you have any questions about your device or its use, please call your JACE Systems Representative.

The following instructions have been prescribed for your CPM use:

TO CHANGE PROGRAM:

1. Stop motion in comfortable position.
2. Press 'PR' arrow twice.
3. Press 'UP' or 'DN' arrow to desired program.
4. Press 'NX' arrow twice. Push GO.

TO CHANGE RANGE OF MOTION:

1. Stop motion in comfortable position.
2. Press 'PR' arrow once.
3. Press 'UP' arrow once.
4. Press 'NX' arrow twice.
5. Extn. Limit will display. Press 'UP' or 'DN' arrow to desired limit.
6. Press 'NX' arrow. Flex. Limit will display. Press 'UP' or 'DN' arrow to desired limit.
7. Press 'NX' arrow (9 times) until screen displays 'Push GO To Start'.
8. Push 'GO'.

NOTE: If the display ever reads 'SET ORTHO UP', press and hold the up arrow until it reads 'Push GO to Start'.