





J • A • C • E  
TRI /// BRACE



**The Integrated CPM Bracing System For Full Knee Extension**



Zero Anterior Tibial Translation ...  
The standard test for AP Translation confirms Zero Anterior Tibial Translation during JACE TRI /// BRACE / K200 operation

Tested by:  
Russell Payne, PT  
Rehab Services of Houston  
Houston, Texas

*"A unique "Dial Locking System" does away with cumbersome pins, screws, hex keys, tool kits or complex set-up procedures. You simply "dial-in" your desired ROM setting."*



Full Knee Extension CPM Brace



Unique Dial Locking System



Post-Operative Brace



Rehabilitation Brace

The **JACE TRI /// BRACE** attaches directly to the **JACE K200 Knee CPM** and is worn during CPM therapy. A convenient interlocking mechanism allows for easy insertion and removal of the leg from the CPM device while maintaining accurate alignment and without the need for additional measurements or realignment when therapy is resumed.

With the **JACE TRI /// BRACE** you can now obtain full knee extension while maintaining the benefits of Zero Anterior Tibial Translation and precise anatomical alignment. The **JACE TRI /// BRACE /K200** eliminates the problem of malalignment often found with other knee CPM devices by ensuring that the knee axis maintains precise alignment with the CPM device axis. Knee stress or torque due to improper alignment is not a problem with **JACE TRI /// BRACE /K200**.

In addition to use as a CPM brace, **JACE TRI /// BRACE** functions as a Post-Operative Brace. **TRI /// BRACE** accommodates both pre-positioned immobilization (0, 30 or 90 degrees) and controlled range of motion requirements (0-120 degrees in 10 degree increments) depending upon your therapy protocol.

The **TRI /// BRACE** easily converts to a Rehab brace. It combines stability and comfort with a full range of protected motion during the latter stages of therapy.





**J · A · C · E**  
**TRI /// BRACE**

The new **JACE TRI /// BRACE** revolutionizes the process of knee rehabilitation therapy with:

- *Full knee extension capability for aggressive rehab protocols*
- *Unique "Dial Locking System"*
- *Direct attachment to the CPM*

The **JACE TRI /// BRACE** is the ideal brace for each stage of your patient's recovery. It functions as a CPM Brace, a Post-Operative Immobilization Brace, and a Rehabilitation Brace. There is no longer a need to incur the expense of additional braces.





Full Knee Extension CPM Brace



Unique Dial Locking System

## JACE TRI//BRACE

- A direct attachment CPM Brace as well as a Post-operative and Rehabilitation Brace
  - Accommodates a variety of patients and therapy needs
  - Full knee extension \*
  - Zero anterior tibial translation \*
  - Precise anatomical alignment \*
  - Custom fit
  - Adjustable, comfortable, safe
  - Ease of application for patient and hospital staff
  - Unique "Dial Locking System" easily adjusts desired ROM setting
  - No pins, screws, hex keys, tools or complex set-up procedures
  - Easy in/out mechanism \*
  - Easily converts to Rehab (ROM) brace
  - Protects against anterior and posterior shift \*
  - Pre-positioned Immobilization at 0°, 30° and 90°
  - Controlled Range of Motion from 0°-120° in 10° increments
- \* With K200 CPM



Post-Operative Brace



Rehabilitation Brace

**J·A·C·E**™  
**SYSTEMS**  
SUPERIOR BY DESIGN

For Additional Information Call  
1-800-800-4CPM  
856-669-6600